

THE CRYOPRO® SERIES OF CRYOSURGICAL TREATMENT UNITS PROVIDES OUTSTANDING PERFORMANCE IN TERMS OF SAFETY, EASE OF OPERATION AND RELIABILITY.

FURTHER, THE CRYOPRO® OFFERS FULL COMPATIBILITY WITH OTHER POPULAR LIQUID NITROGEN UNITS, I.E. SPRAY APERTURES AND CONTACT PROBES FOR OTHER SYSTEMS MAY BE USED WITH THE CRYOPRO® UNITS AS WELL.

THE LOW TEMPERATURE, THE VARIETY OF ACCESSORIES FOR BOTH SPRAY AND CONTACT FREEZING AND THE HIGH CURE RATE COMBINED WITH EXCELLENT COSMETIC RESULTS MAKES LIQUID NITROGEN CRYOSURGERY THE MOST VERSATILE CRYOSURGICAL TREATMENT AVAILABLE AT A VERY LOW COST PER TREATMENT.

THE CRYOPRO® IS AVAILABLE IN TWO SIZES:

CRYOPRO® MAXI WITH 500 ML CAPACITY (20 - 24 HRS. STATIC HOLDING TIME) AND **CRYOPRO® MINI** OFFERING A 350 ML CAPACITY (12 - 14 HRS. STATIC HOLDING TIME). FURTHER, A COMPLETE RANGE OF ACCESSORIES IS AVAILABLE INCL. SPARY TIPS, CONTACT PROBES, DEWARS, WITHDRAWAL DEVICES ETC. THREE YEAR WARRANTY.



CORTEX TECHNOLOGY



Smedevaenget 10 9560 Hadsund Denmark
Tel.: +45 9857 4100 Fax: +4598572223
E-mail: cortex@cortex.dk WEB: www.cortex.dk

German office:
Tel.: +49 2572 941046 Fax: +49 2572 941047
E-mail: cortexde@aol.com



CryoPro®

Liquid Nitrogen Cryosurgery



“CryoPro® is not cheap, but it has brought some significant improvements to our cryotherapy usage.”

■ ■ ■

“The CryoPro® is infinitely better than the hit-or-miss approach of a cotton wool bud.”

Dr. G. Caranci (GP Medicine, January 15 1999)



“In routine clinical practice, liquid nitrogen spray equipment has become increasingly dominant over other methods.”

Dawber, Colver, Jackson: Cut.
Cryosurgery (Martin Dunitz)